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GB INSTRUCTIONS FOR USE BIOTECH DENTAL DENTAL IMPLANT SYSTEMS

KONTACT®, B.I.S.®, B.I.S. Conic® and Weego.fr.

These products must be handled by TRAINED, QUALIFIED individuals who HAVE READ these instructions.

1. GENERAL DESCRIPTION

BIOTECH DENTAL dental implant systems consist of:

- titanium/titanium alloy dental implants which comply with applicable standards ISO 5832-2/ISO 5832-3,
- titanium, titanium alloy and ceramic abutments which comply with the applicable standards,
- additional prosthetic titanium components in titanium, titanium alloy or metal and polymer which comply with the applicable standards,
- ancillary medical grade stainless steel instrumentation, titanium/titanium alloy. For the purpose of identification, the majority of the instruments bear a visible reference and manufacturing batch number. The sterile products are identified on the packaging. If the word 'sterile' is not specifically stated, the products are supplied non-sterile and must be sterilised before use by the staff concerned. The dental implant systems are radio-opaque.

2. GENERAL CONDITIONS OF USE

- These products must only be used by experienced surgeons;
 - trained in implantology and techniques for the use of BIOTECH DENTAL medical devices,
 - in the observance of the indications, contraindications and precautions which are to be indicated to the patient during the postoperative phase and which are mentioned in these instructions as well as in surgical and prosthetic protocols available on request,
 - on premises compatible with the aseptic conditions required by these instructions.
- Under no circumstances may the dental implant system be used with the devices of another manufacturer. During the period from implantation to management, BIOTECH DENTAL undertakes to replace any implant or prosthetic component, provided that the above conditions have been observed. The practitioner will be responsible for complications which result from use which has failed to comply with the recommendations or asepsis. These complications may under no circumstances be attributed to BIOTECH DENTAL and will lead to the rejection of any formulated action against BIOTECH DENTAL concerning guarantee or product replacement.

3. INDICATIONS

BIOTECH DENTAL dental implants are used as a base for the anchorage of replacements to the mandible as well as the maxillary. Restorations range from unitary replacement to complete replacement for the purposes of reproducing all or some of the functions of the natural teeth. This technique is only considered if there less invasive therapies are regarded as less suitable. The osseous volume and quality must be sufficient for the purposes of supporting the dental implants.

4. CONTRAINDICATIONS

Dental implants must not be used in patients who present with local and/or general contraindications to surgical dental intervention. These contraindications must be carefully assessed when the surgeon makes his decision. Among other things, the contraindications constitute:

- factors which interfere with the process of healing and osseous and tissue regeneration such as problems with osseous metabolism, connective tissue, type 2 diabetes, steroid therapy, osseous infections, and nicotine and alcohol abuse,
- diseases or immunosuppressive therapies such as chemotherapy and radiotherapy,
- osseous volume and osseous quality which are not sufficient for the purposes of accepting a sufficient number and size of implants and to support the functional loads,
- bruxism, parafunctional habits, problems with occlusion and/or the temporomandibular joints,
- buccal infections and inflammation such as parodontitis and gingivitis,
- allergy to titanium or titanium alloy or inadequate bucco-dental hygiene.

5. PRECAUTIONS FOR USE

- In order to avoid swallowing or inhalation of small components, it is recommended that they be secured where possible.

- BIOTECH DENTAL must be informed of any incident or problem concerning an implant or instrument.

Preoperative precautions :

- The user must ensure that the material is in good condition and working correctly before use. Under no circumstances must the material be used if the packaging is damaged or if surface damage or deformations can be seen.
- Evaluation of the patient is necessary for the purposes of determining any factor which is likely to expose him to a risk resulting from insertion of the implant or which might have an influence on healing of the bone/soft tissue. Certain criteria must be taken into consideration with regard to the indication for surgery and the choice of implant.

- The occlusal surface of the implant must ideally be smaller than the prosthetic tooth to ensure the flaring of the soft tissue and the emergence profile of the prosthesis.

- a minimum space of 1.5 mm between the implant and the adjacent root, 1.5 mm between 2 implants, and a thickness of bone of 1 mm around the vestibular and lingual surfaces of the implant (> in the case of large diameters),

- in the event of suspected hypersensitivity or proven hypersensitivity to foreign bodies, it is recommended that titanium and TARGV ELI titanium tolerance be checked prior to implantation of the material.

Radiographic and/or CT scan evaluation and a clinical treatment plan are essential for the purposes of ensuring complete safety.

Postoperative precautions :

- Primary stability of the implant is essential,
- The choice of the type of implant and its size will be a matter for the surgeon; this will depend on the particular patient and the patient's pathology. The size and the nature of the osseous structures will determine the surgeon's choice. The surgeon will also have to consider the stresses on the implant after the procedure,
- Implants must not be deformed. At the risk of causing abnormal fatigue of the implant, the manipulations will have to be progressive. An excessive effort exerted on the implant might induce stresses which could cause a fracture or deformation of the implant responsible for secondary effects,
- It is up to the practitioner to determine the settings of his equipment depending on the particular clinical case, and to check that his equipment is working properly prior to every procedure. It is recommended that the drilling sequences be conducted at a slow speed and with profuse and constant external irrigation. The recommended rotation speed is from 700 to 1,200 rpm in the case of drilling, not exceeding 1,500 rpm, and 200 to 400 rpm for boring,
- The drills are very sharp; it is recommended that caution be exercised when gripping the drills in order to avoid being cut. Also, because the drills are designed with helicoidal flutes, their use can result in a drag phenomenon which is detrimental to drilling precision. In order to correct this possibility, it is recommended that use be made of the drilling stops provided for this purpose.
- Cutting tools (countersinks, drills, taps and borer)s must not be used more than 20 times in order to prevent problems with wear. A specific tracking sheet for the drills is available.

Postoperative precautions :

After insertion of the implant, osteointegration is achieved within 4 months (mandible) 6 months (maxillary). The surgeon must implement postoperative follow-up. He will be responsible for informing the patient with regard to actions to avoid following the insertion of an implant and the precautions to be taken in the course of normal daily activity. The patient must understand that prior to complete osseous consolidation, a metallic material is not as strong as healthy osseous structures, and that metallic materials will break if subjected to intense or abnormal stresses.

- the efficacy of the implant depends on his behaviour and his buccal hygiene,
- potential risks and adverse effects associated with the insertion of a dental implant,
- the possibility of failure of implant osteointegration,
- the need to report any sensitivity appearing in the operated area, pre- or post-operatively,
- potential interaction with other procedures and therapeutic or diagnostic devices.

6. POSSIBLE ADVERSE EFFECTS

The implantation techniques are accompanied by contraindications and risks of use. They can also be related to use of the device as well as the surgical intervention and they can include (in a non-exhaustive list):

- loosening, rupture or loss of a component,
- local nervous lesion, transient or permanent,
- infection, lesion of adjacent teeth and aesthetic problems,
- an allergic reaction to the different materials which constitute BIOTECH DENTAL implants.

7. STERILITY AND TRACABILITY

Single delivered device : dental implants delivered with a cover screw and a healing screw in the packaging of the implant, are sterile (sterilisation by gamma radiation). Any product which is supplied sterile is for single use and must not be re-sterilised. The entire packaging must be checked before use to ensure that the sterility of the contents has not altered. In conditions of adequate storage, the packaging maintains sterility until the expiry date. Do not use if the packaging has been damaged.

Labels on the different levels of packaging as well as free labels will make it possible to ensure perfect traceability based on the reference number and the batch number of the implant.

Single-use single-use devices should under any circumstances be reused and re-sterilized. In case of re-sterilization of a non-re-sterilizable medical device, the absence of iatrogenic transmission of pathogens by these devices is no longer guaranteed and accident in the field of materiogilance may occur.

Non-sterile delivered device : the ancillary instrumentation and the prosthetic parts are supplied clean but not sterile. It will be the responsibility of the staff concerned to proceed with decontamination, cleaning and sterilisation before and after every procedure (for reusable products) in accordance with the applicable recommendations.

8. CLEANING AND STERILISATION OF NON-STERILE DEVICES

Decontamination and Cleaning operations are imperative before every sterilisation process.

- In order to prevent the deterioration of the components, only use cleaning products and decontamination products which are compatible with the materials. Alkaline and corrosive products for stainless steel and aluminium medical devices are prohibited; bleach, aqueous solutions of sodium hypochlorite and sodium chloride (bleach) and/or corrode the instruments,
- the water to be used for pre-désinfection, cleaning, decontamination, rinsing and sterilisation must comply with the applicable regulations,
- The user must pay particular attention to the cleaning of dirt, residue and deposits on all parts of the instruments (cavities, interstices, etc.).

Pre-decontamination :

This must be done immediately after every procedure in the case of all disassembled reusable instruments, completely immerse the previously disassembled instruments in a decontaminating solution with a neutral pH, without aldehyde and whose concentration depends on the rate of contamination and the degree of soiling of the instruments. Soak for 15 to 20 minutes. Ultrasonic bath soaking is recommended,

- the decontamination solution must be renewed each time after use in order to prevent its saturation.

Cleaning :

- This can be undertaken by brushing and by ultrasound and must be done on disassembled components,
- brushing must be carried out with a soft brush (for example, nylon),
- soaking for a minimum of 15 minutes is necessary. Prolonged contact with certain products for several hours can cause corrosion of the material (do not exceed 2 hours).

Rinsing – Drying :

Cleaning must immediately be followed by abundant rinsing and then the products must be dried prior to beginning the sterilisation process.

Packaging prior to sterilisation :

Whichever method of cleaning is used, the material must be packaged immediately in order to avoid any new contamination. In the event of sterilisation on metal grid trays, double wrapping paper must be used for the packaging. In the event of late packaging (> 2 hours), further cleaning might be considered in order to make sure that the material is sterile.

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Prepare the instruments so that all surfaces are in direct contact with the steam. The recommended cycle is the P101 cycle at 134°C for 18 minutes.

Retour :

After using the ancillary in the event of a loan, this must be returned to BIOTECH DENTAL after having been cleaned and decontaminated in accordance with applicable regulations. It is incumbent upon the healthcare establishment or the practitioner to ensure that the cleaning/sterilisation process which is effectively implemented (equipment, materials and operator) achieves the desired result. Under no circumstances will the manufacturer be held responsible for sterility.

9. HANDLING AND STORAGE

The dental implant system must be handled and stored carefully in order to ensure a lack of contamination and deterioration prior to any surgical intervention.

Storage must be done with care in an appropriate, dry and clean environment. They must not be exposed to direct sunlight, ionising radiation, extreme temperatures, specific contamination, or contact or proximity with products that may have a corrosive effect.

The lifecycle of the devices can be compromised by a lack of precautions during handling or poor protection of the instruments.

10. ELIMINATION

The elimination of medical waste must comply with the applicable legislation. The professional who produces waste will be responsible for sorting waste which carries a risk of infection and waste which is compatible with household waste.

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